

UCSF Medical Center
Point of Care Testing

Clinitek Status Maintenance LOG

Location: _____ Month/Year: _____ Serial Number # _____

Date	Wipe outside with damp cloth/mild detergent. Wipe Test table insert. <u>Daily</u>	Test Table insert cleaning. Inspect calibrator for scratches <u>Weekly</u>-(first work day each week)	Disinfecting the test table and insert. <u>Monthly</u>-(first week of the month)	Comments/ Corrective Action	Initials
1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
11					
12					
13					
14					
15					
16					
17					
18					
19					
20					
21					
22					
23					
24					
25					
26					
27					
28					
29					
30					
31					

Reviewed By: _____ Date: _____