

SAMPLE COLLECTION GUIDE FOR ITL TESTS:

Please use the following guide when collecting sample for ITL:

ITL TEST CATEGORY	MAXIMUM # TUBES FOR PT's > 12 KG	SAMPLE VOLUME FOR PT's < 12 KG
Antibody	2 RED TOP	3 mL in RED TOP
HLA Typing	2 ACD*	6 mL in small ACD
Chimerism	2 ACD*	6 mL in small ACD
Crossmatch Donor	6 ACD*	CONTACT ITL
Crossmatch Recipient	1 RED TOP	CONTACT ITL
Auto Crossmatch	6 ACD* & 1 RED TOP	CONTACT ITL

* Completely full 8.5 mL tubes

If more than one test is ordered in the same category, collect no more than the maximum tube or sample volume indicated on the chart.

If two or more test categories are ordered, add the # tubes together.

ITL (415) 476-3886