

Machaon Diagnostics

ADAMTS13 Test Order Form

National Service with Labs in California and Louisiana
 Medical Director: Brad H. Lewis, MD (Oakland, CA lab)
 Medical Director: Gloria Coker, MD (New Orleans, LA lab)
 Phone: (510) 839-5600 Fax: (510) 839-6153

TAT <24 hours, 7 days a week; order STAT for weekend testing

MachaonDiagnostics.com
 1-800-566-3462

PATIENT INFORMATION (complete or attach)			SUBMITTING FACILITY		
Patient's Name: (Last, First, M.I.)		Gender:	UCSF Clinical Labs 505 Parnassus Ave., Room M524 San Francisco, CA 94143 Verify receipt of ALL faxed results: (415) 353-1747 PARNASSUS HEMATOLOGY LAB. Fax results to China Basin (415) 353-1804 and Parnassus (415) 514-5826. Call all specimen problems to China Basin Lab at (415) 353-1349.		
Specimen Date: (MM/DD/YY) Time: AM / PM		Date of Birth: (MM/DD/YY)			
Medical Record / ID#:		Patient's Social Security #:			
ORDERING PHYSICIAN INFORMATION			BILLING INFORMATION		
Physician's Name: (Last, First, M.I.)			Bill to: <input type="checkbox"/> Facility <input type="checkbox"/> Insurance <input type="checkbox"/> Patient <input type="checkbox"/> Medicare		
Physician's direct phone number to call results: (called within 24 hours)			Patient status: <input type="checkbox"/> Inpatient <input type="checkbox"/> Outpatient		
Fax Number for Results:			<input type="checkbox"/> STAT <i>Samples shipped for weekend analysis must be ordered STAT. Mark FedEx Airbill for 'SATURDAY Delivery.'</i>		
CLINICAL INFORMATION (if available)			TEST SELECTION		
Last Plasma Infusion Date: (MM/DD/YY)		LDH: (U/L)	<input type="checkbox"/> ADAMTS13 Activity reflex Inhibitor		<i>Reflexes to Antibody</i>
Creatinine: (mg/dL)	PLT Count: (K/ μ L)	Hemoglobin: (g/dL)	<input type="checkbox"/> ADAMTS13 Activity, Inhibitor and Antibody		
Clinical Suspicion: (circle one) TMA / TTP / HUS / Shiga toxin-related HUS / aHUS / Other (Call lab with ADAMTS13 Gene Sequencing requests)			<input type="checkbox"/> Reflex to aHUS Genetic Panel when ADAMTS13 Activity is >5% (3mL EDTA whole blood sample required)		
SPECIMEN COLLECTION, PROCESSING AND SHIPPING / call for draw kits (free shipping)					
Specimen Collection and Processing			Specimen Shipping		
1. Draw 1 blue top tube (3.2% sodium citrate).			6. Place all forms into document sleeve of the biohazard bag.		
2. Spin tube for cell-free plasma (approx. 15 minutes at 2000 x g).			7. Ensure that no patient-specific information is visible.		
3. Transfer approximately 1mL into the provided sample tube.			8. Place biohazard bag into the provided FedEx Clinical Pak.		
4. Label tube with patient first and last name, draw date and DOB.			9. Complete the provided FedEx Airbill and affix Airbill pouch to Pak.		
5. Enclose tube in Styrofoam box and seal in biohazard bag.			10. Call FedEx 800-238-5355 for a pick-up.		
Note: Samples can be shipped room temperature, frozen or refrigerated. Please call for guidance.					
OUTPATIENT ONLY: INSURANCE BILLING INFORMATION (complete or attach)					
Insurance Company: (Medicare patients must sign ABN on reverse)			Patient Address:		Patient Phone Number:
Insurance Policy / Medicare Number:		Insurance Group Number:	Patient City:		State: Zip Code:
Insurance Company Address:			Authorization Number:		
Insurance Company City:		State:	Zip Code:		
DIAGNOSIS CODE(S):			OUTPATIENT ONLY: PATIENT SIGNATURE		
ICD-10 Code:			Machaon Diagnostics may need to obtain additional information from your physician to complete these services. I hereby authorize the release of medical information related to the services described herein and authorize payment directly to Machaon Diagnostics. Machaon Diagnostics is a PARTICIPATING PROVIDER of Medicare only. The activity test is \$248 and reflexes to the inhibitor (\$530) and antibody test (\$250); shipping charges may apply. I agree to assume responsibility for payment of all charges not covered by my healthcare insurer.		
ICD-10 Code:					
ICD-10 Code:			Patient's Signature: _____ Date: _____		
ADDITIONAL INFORMATION					

Machaon Diagnostics is a specialized coagulation and platelet laboratory that provides comprehensive clinical evaluations of bleeding and clotting patients. Most evaluations can be completed within 24 hours, 7 days a week. Machaon Diagnostics is a California-licensed, CLIA-accredited, CAP-accredited, clinical laboratory approved to provide high-complexity testing services. Machaon Diagnostics is a PARTICIPATING PROVIDER of Medicare only. Patients with insurance coverage other than Medicare are considered OUT-OF-NETWORK and will be billed for services not covered by their insurance provider. Medicare patients must sign an ABN, either located on the reverse side of this form or downloaded from the Machaon Diagnostics website. Patient insurance billing services are provided in accordance with the Machaon Insurance Billing Policy. HMO or medical group covered patients may need a prior authorization if they seek full reimbursement. For more information please visit www.MachaonDiagnostics.com or call (510) 839-5600.

MDI Use: (Order number): _____ (Number of aliquots): _____

Version:19OCT2020